

mainwindow

```
+ bool flag

- UI::MainWindow *ui
- Device* device

- int selectedSession
- int selectedGroup
- int connectionIntensity
- int timesIntensityAdjusted
- bool endSessionEarlyFlag

- QElapsedTimer blinkTimer
- QElapsedTimer powerPressedTimer
- QElapsedTimer selectTimer
- QElapsedTimer pauseTimer
- QElapsedTimer therapyTimer
- QElapsedTimer timeoutTimer
- QElapsedTimer testingTimer

- therapy(int groupNum, int sessionNum, int recordingFlag, int overrideIntensity)
- connectionTest()
- addRecording(const string& name, int group, int batteryPercent, int initialIntensity, int intensity)
- replayRecording(Recording* recording)


- changeTextColor(QTextBrowser *text, QColor color)
- changeBackgroundColor(QPushButton *button, const QString& color, const QString& image, const QString& radius)
- bootConnectionTest()
- blinkSession(int sessionNum)
- cycleGroupButton()
- turnOffUI()
- softOff()
- setTherpayLock(bool status)
- updateIntensityLog()
- displayIntensityOnGraph()
- powerLightOnOff(bool status)
- endSessionEarly()

- ui_initializeBattery()
- ui_turnOffBattery()


- blinkBattery()
- drainBattery(int intensity)
- checkBattery(): bool
- batteryWarning()

- sleepy(int sleepTime)
```






Power On




User presses the power button starting the device. Device ensures the battery is at acceptable levels, if it is the device starts and lights are turned on. If it isn't the user is prompted to change batteries. Refers to the Power On use case (UC1).

Power Off



The user presses and holds the power button, signaling the device to shutoff. If it is running a session the intensity is slowly lowered before shutting off. Alternatively, if the device has been idle for too long after booting the device shuts off. Refers to the Power Off use case (UC2).

Low Battery




The main window periodically checks the battery level.

If battery levels are sufficient functioning continues, else the battery light blinks and user is prompted to change batteries.

Refers to the Low Battery use case (UC3)

Selecting a Session




User wants to select a session, presses the up/down buttons to change selection.

Main window is updated as to which selection is active, lighting up on the device.

After user confirms session, it alerts the user then starts after 5 seconds. Refers to the Selecting a Session use case(UC4)


Connection Test



The connection test begins upon startup or if the device is disconnected during a session. If the device is disconnected, buttons will be disabled until a connection is re-established. The mid-session alternative shows the session being frozen until a connection is re-established. If not mid-session, the device blinks the intensity lights to indicate it is looking for a connection. It continues this behaviour until a connection is established.

This refers to the connection test use case (UC5)

Intensity




When the user presses the up/down arrow button during or while configuring a session, the intensity is increased or decreased.

The device first checks that the desired intensity is within the range or possible intensities and if not, a warning is displayed.

This refers to the Intensity use case (UC6)


Record



If the user holds the select button for a second they prime the device to remember the current settings. When this happens, the user can choose to press 'Stop' which saves the settings but does not execute the session. Otherwise the session's settings are saved and the session is executed.

This refers to the Record use case (UC7)

Replay



If the user presses the Print History Button, the device goes through the user's recording sessions and displays their settings.

If the user has no recorded sessions the device displays "No History".

Otherwise the user may select a recorded session, then press select to begin a session with the recorded settings.

This refers to the Replay use case (UC8)